

## Implant Summary Report for:

Smith & Nephew

### R3 Cementless Cup

Comprising PRIMARY hips implanted up to:  
NJR Database extract:

09 September 2018  
08 November 2018

Produced on:  
Licenced for use until:

18 November 2018  
18 March 2019

## Contents

Recorded Usage in NJR	2
Patient and Procedure Details	3
Revision and Survivorship	4
Patient Reported Outcomes	12
APPENDIX A – Component List	15

This report has been produced by the National Joint Registry of England, Wales and Northern Ireland (NJR). It summarises usage and outcomes associated with the R3 Cementless Cup up to the specified dates. This analysis is based on data collected by the NJR and PROMs data collected by NHS Digital for the following combinations of products:

**Acetabular components**

**R3 Cementless Cup**

**Femoral Components**

**Any Femoral Components manufactured by Smith & Nephew**

Appendix A lists details of all R3 Cementless Cup components included in this report.

This report is made available to Smith & Nephew under licence. Under this licence, Smith & Nephew is free to distribute this document, subject to the restrictions specified in the National Joint Registry Supplier Feedback Data Access and Use Code of Practice.

**NJR Recorded Usage**

**Implant Usage**



Date of first recorded usage in the NJR: **27 March 2007**

Date of last recorded usage in this dataset: **08 September 2018**

Maximum implantation time: 11.4 years

Mean implantation time: 3.2 years

Totals Recorded in NJR	Cumulative Total
Procedures	22,030
Patients	20,053
Centres	126
Consultants	312
Implanting Surgeons	504



Current Outcome	Year of implantation														Total	%
	Pre-2007	2007	2008	2009	2010	2011	2012	2013	2014	2015	2016	2017	2018			
<b>Death</b>	0	1	10	38	63	80	110	185	147	86	67	39	7	833	3.8%	
<b>Revised</b>	0	5	13	16	11	21	31	35	43	55	24	31	11	296	1.3%	
<b>Unrevised</b>	0	6	52	277	511	872	1530	2282	2638	2893	3216	3682	2942	20901	94.9%	
<b>Total</b>	<b>0</b>	<b>12</b>	<b>75</b>	<b>331</b>	<b>585</b>	<b>973</b>	<b>1671</b>	<b>2502</b>	<b>2828</b>	<b>3034</b>	<b>3307</b>	<b>3752</b>	<b>2960</b>	<b>22030</b>	<b>100%</b>	

## Patient / Procedure Details

## Patient Details

	R3 Cementless Cup	All NJR Cementless Cups	All NJR hip replacement
<b>Total Procedures</b>	22,030	678,283	1,070,534
<b>Total Patients</b>	20,053	582,437	912,154
<b>Demographics</b>			
Mean age	<b>65.0</b>	<b>65.9</b>	<b>68.6</b>
< 50	10.1%	8.5%	6.0%
50 – 59	19.5%	18.2%	13.7%
60 – 69	32.3%	32.7%	29.2%
70 – 79	28.9%	29.3%	34.6%
≥ 80	9.2%	11.2%	16.4%
Median BMI	<b>28</b>	<b>28</b>	<b>28</b>
% BMI information available	74.8%	61.8%	59.7%
Underweight ( BMI < 18.5)	0.6%	0.8%	0.9%
Normal (18.5 ≤ BMI < 25)	18.1%	19.9%	20.8%
Overweight (25 ≤ BMI < 30)	39.3%	39.5%	39.7%
Obese I (30 ≤ BMI < 35)	27.3%	26.1%	25.7%
Obese II (35 ≤ BMI < 40)	11.0%	10.1%	9.6%
Obese III (BMI ≥ 40)	3.6%	3.5%	3.3%
% Male	44.5%	43.8%	40.2%
<b>ASA Grades</b>			
P1 - Fit and healthy	18.0%	19.1%	16.3%
P2 - Mild disease not incapacitating	67.6%	67.3%	67.6%
P3 - Incapacitating systemic disease	14.1%	13.2%	15.6%
P4 / P5	0.3%	0.4%	0.6%
<b>Indications</b>			
Osteoarthritis	91.88%	91.97%	91.72%
Rheumatoid Arthritis	1.43%	1.28%	1.32%
Avascular Necrosis	2.46%	2.66%	2.51%
Fractured Neck of Femur	2.39%	2.72%	3.20%
CDH/DDH	2.66%	2.07%	1.55%
Other	1.55%	1.99%	2.13%

Total of indications may exceed total number of implants, since more than one indication can be listed per case.

## Revision and Survivorship

## Reasons for Revision - All bearing types

Reasons for revision of primary procedures in which the implant was used.

Reason for Revision	Revised <sup>†</sup>	Expected Revisions <sup>*</sup>	p value
Unexplained Pain	26	74.77	<0.001
Dislocation / Subluxation	69	77.80	0.358
Adverse Soft Tissue Reaction	24	53.06	<0.001
Infection	57	65.45	0.317
Aseptic Loosening - Stem	39	55.47	0.025
Aseptic Loosening - Socket	10	38.95	<0.001
Periprosthetic Fracture Stem	49	57.38	0.286
Periprosthetic Fracture Socket	3	7.77	0.099
Malalignment Stem	10	13.40	0.408
Malalignment Socket	14	26.38	0.013
Wear Of Acetabular Component	13	13.35	1
Lysis Stem	3	8.90	0.042
Lysis Socket	1	7.89	0.006
Implant Fracture Stem	6	6.08	1
Implant Fracture Socket	1	6.78	0.017
Implant Fracture Head	1	2.80	0.536
Dissociation of Liner	3	7.65	0.098
Other / reason not recorded	17	29.23	0.019
<b>Total Revised</b>	<b>296</b>	<b>433.44</b>	<b>&lt;0.001</b>

† multiple reasons may be listed for one revision procedure

\* Based on All NJR Cementless Cups adjusted for agegroup, gender, and indications

Significantly better,  $p < 0.001$ Significantly better,  $p < 0.05$ Significantly worse  $p < 0.05$ Significantly worse  $p < 0.001$ 

Components Revised	Number of procedures	R3 Cementless Cup	All other Cementless Cups in NJR
Femoral only	99	33%	21%
Acetabular only	47	16%	26%
Both femoral and acetabular	79	27%	44%
Neither femoral nor acetabular revision recorded *	71	24%	9%

\* includes isolated head and/or liner exchange

## Revision and Survivorship

## Reasons for Revision - excluding Metal on Metal

Reasons for revision of non-MoM primary procedures in which the implant was used.

Reason for Revision	Revised <sup>†</sup>	Expected Revisions <sup>*</sup>	p value
Unexplained Pain	21	40.22	0.001
Dislocation / Subluxation	66	78.53	0.169
Adverse Soft Tissue Reaction	6	9.31	0.404
Infection	53	63.81	0.182
Aseptic Loosening - Stem	38	48.90	0.128
Aseptic Loosening - Socket	9	25.16	<0.001
Periprosthetic Fracture Stem	48	50.66	0.776
Periprosthetic Fracture Socket	3	7.30	0.13
Malalignment Stem	10	12.38	0.663
Malalignment Socket	12	22.16	0.03
Wear Of Acetabular Component	12	11.71	0.882
Lysis Stem	3	5.11	0.501
Lysis Socket	1	3.77	0.195
Implant Fracture Stem	6	5.20	0.655
Implant Fracture Socket	0	7.69	<0.001
Implant Fracture Head	1	2.36	0.735
Dissociation of Liner	3	7.86	0.099
Other / reason not recorded	13	18.05	0.281
<b>Total Revised</b>	<b>263</b>	<b>338.48</b>	<b>&lt;0.001</b>

† multiple reasons may be listed for one revision procedure

\* Based on All NJR Cementless Cups (excluding metal on metal), adjusted for agegroup, gender, and indications


Significantly better,  $p < 0.001$ Significantly better,  $p < 0.05$ Significantly worse  $p < 0.05$ Significantly worse  $p < 0.001$ 

Components Revised	Number of procedures	R3 Cementless Cup	All other Cementless Cups in NJR
Femoral only	98	37%	30%
Acetabular only	40	15%	25%
Both femoral and acetabular	69	26%	34%
Neither femoral nor acetabular revision recorded *	56	21%	11%

**Revision and Survivorship**

**Cumulative Revision Rate**

Endpoint: All reasons for revision. All Bearing types




Cox Proportional Hazards model for revision risk ratio of R3 Cementless Cup / All other Cementless Cups in NJR, with endpoint as any revision.

Adjustment	Hazard Ratio (95% CI)	p-value
All bearings. Unadjusted.	0.69 (0.62 - 0.77)	<0.001
All bearings. Adjusted for age, gender, year cohort and indications.	0.88 (0.78 - 0.99)	0.032

**Revision and Survivorship**

**Cumulative Revision Rate**

Endpoint: All reasons for revision. Excluding Metal-on-Metal




Cox Proportional Hazards model for revision risk ratio of R3 Cementless Cup / All other Cementless Cups in NJR, with endpoint as any revision.

Adjustment	Hazard Ratio (95% CI)	p-value
Excluding MoM, Unadjusted.	0.79 (0.70 - 0.89)	<0.001
Excluding MoM. Adjusted for age, gender, year cohort and indications.	0.83 (0.73 - 0.94)	0.003

**Revision and Survivorship**

**Cumulative Revision Rate**

Endpoint: Acetabular Revision. All Bearing types



Cox Proportional Hazards model for revision risk ratio of R3 Cementless Cup / all other All other Cementless Cups in NJR, with endpoint as Acetabular revision.


Adjustment	Hazard Ratio (95% CI)	p-value
All bearings. Unadjusted.	0.45 (0.38 - 0.54)	<0.001
All bearings. Adjusted for age, gender, year cohort and indications.	0.65 (0.55 - 0.78)	<0.001



**Revision and Survivorship**

**Cumulative Revision Rate**

Endpoint: Acetabular Revision. Excluding Metal-on-Metal



Cox Proportional Hazards model for revision risk ratio of R3 Cementless Cup / All other Cementless Cups in NJR, with endpoint as Acetabular revision.

Adjustment	Hazard Ratio (95% CI)	p-value
Excluding MoM, Unadjusted	0.57 (0.47 - 0.68)	<0.001
Excluding MoM. Adjusted for age, gender, year cohort and indications.	0.63 (0.52 - 0.76)	<0.001

Revision and Survivorship

Cumulative Revision Rate

Endpoint: All revisions

Period /years	All Bearing types				Excluding Metal on Metal Bearings			
	At Risk	R3 Cementless Cup	All NJR Cementless Cups	All THR	At Risk	R3 Cementless Cup	All NJR Cementless Cups	All THR
0	22,030	-	-	-	21,891	-	-	-
1	17,469	0.8% (0.7% - 0.9%)	0.9% (0.9% - 0.9%)	0.8% (0.8% - 0.8%)	17,336	0.8% (0.7% - 0.9%)	0.9% (0.9% - 0.9%)	0.8% (0.7% - 0.8%)
2	13,805	1.1% (0.9% - 1.2%)	1.3% (1.3% - 1.4%)	1.2% (1.1% - 1.2%)	13,677	1.0% (0.9% - 1.2%)	1.2% (1.2% - 1.3%)	1.1% (1.1% - 1.1%)
3	10,435	1.3% (1.1% - 1.4%)	1.7% (1.7% - 1.8%)	1.5% (1.5% - 1.5%)	10,311	1.2% (1.1% - 1.4%)	1.5% (1.5% - 1.6%)	1.4% (1.4% - 1.4%)
4	7,614	1.4% (1.2% - 1.6%)	2.2% (2.1% - 2.2%)	1.9% (1.9% - 1.9%)	7,497	1.3% (1.1% - 1.5%)	1.8% (1.8% - 1.8%)	1.6% (1.6% - 1.6%)
5	4,855	1.7% (1.5% - 1.9%)	2.7% (2.6% - 2.7%)	2.3% (2.3% - 2.3%)	4,742	1.6% (1.4% - 1.8%)	2.0% (2.0% - 2.1%)	1.9% (1.8% - 1.9%)
6	2,780	2.0% (1.7% - 2.3%)	3.3% (3.2% - 3.3%)	2.7% (2.7% - 2.8%)	2,672	1.7% (1.5% - 2.0%)	2.3% (2.3% - 2.3%)	2.1% (2.1% - 2.2%)
7	1,412	2.5% (2.1% - 3.0%)	3.9% (3.8% - 3.9%)	3.2% (3.2% - 3.3%)	1,315	2.0% (1.7% - 2.4%)	2.6% (2.5% - 2.6%)	2.4% (2.3% - 2.4%)
8	677	3.0% (2.5% - 3.6%)	4.5% (4.5% - 4.6%)	3.7% (3.7% - 3.8%)	591	2.2% (1.8% - 2.8%)	2.9% (2.8% - 2.9%)	2.7% (2.6% - 2.7%)
9	226	3.5% (2.7% - 4.5%)	5.2% (5.1% - 5.3%)	4.3% (4.2% - 4.3%)	177	2.2% (1.8% - 3.2%)	3.2% (3.1% - 3.3%)	3.0% (3.0% - 3.1%)
10		(-)	5.9% (5.8% - 6.0%)	4.8% (4.8% - 4.9%)		(-)	3.6% (3.5% - 3.7%)	3.4% (3.3% - 3.4%)
11		(-)	6.6% (6.5% - 6.7%)	5.4% (5.3% - 5.5%)		(-)	4.0% (3.9% - 4.1%)	3.8% (3.8% - 3.9%)
12		(-)	7.3% (7.1% - 7.4%)	6.0% (5.9% - 6.1%)		(-)	4.4% (4.3% - 4.6%)	4.3% (4.2% - 4.4%)
13		(-)	7.9% (7.7% - 8.0%)	6.6% (6.5% - 6.7%)		(-)	4.9% (4.8% - 5.1%)	4.8% (4.7% - 4.9%)
14		(-)	8.6% (8.4% - 8.8%)	7.2% (7.0% - 7.3%)		(-)	5.5% (5.2% - 5.7%)	5.2% (5.1% - 5.4%)
15		(-)	9.4% (9.1% - 9.7%)	7.9% (7.6% - 8.1%)		(-)	(-)	6.0% (5.7% - 6.2%)

Cumulative revision rate with 95% confidence intervals

Rate is only reported for times where more than 40 remain at risk

Revision and Survivorship

Cumulative Revision Rate

Endpoint: Acetabular Revision

Period /years	All Bearing types				Excluding Metal on Metal Bearings			
	At Risk	R3 Cementless Cup	All NJR Cementless Cups	All THR	At Risk	R3 Cementless Cup	All NJR Cementless Cups	All THR
0	22,030	-	-	-	21,891	-	-	-
1	17,469	0.3% (0.3% - 0.4%)	0.5% (0.5% - 0.5%)	0.5% (0.4% - 0.5%)	17,336	0.3% (0.2% - 0.4%)	0.5% (0.5% - 0.5%)	0.4% (0.4% - 0.4%)
2	13,805	0.4% (0.3% - 0.5%)	0.8% (0.8% - 0.8%)	0.7% (0.7% - 0.7%)	13,677	0.4% (0.3% - 0.5%)	0.7% (0.7% - 0.7%)	0.7% (0.6% - 0.7%)
3	10,435	0.5% (0.4% - 0.6%)	1.1% (1.1% - 1.1%)	1.0% (1.0% - 1.0%)	10,311	0.5% (0.4% - 0.6%)	0.9% (0.9% - 0.9%)	0.9% (0.8% - 0.9%)
4	7,614	0.6% (0.5% - 0.7%)	1.4% (1.4% - 1.4%)	1.3% (1.2% - 1.3%)	7,497	0.5% (0.4% - 0.7%)	1.0% (1.0% - 1.1%)	1.0% (1.0% - 1.0%)
5	4,855	0.8% (0.6% - 0.9%)	1.8% (1.8% - 1.8%)	1.6% (1.5% - 1.6%)	4,742	0.7% (0.5% - 0.8%)	1.2% (1.2% - 1.2%)	1.2% (1.2% - 1.2%)
6	2,780	0.9% (0.7% - 1.1%)	2.3% (2.2% - 2.3%)	1.9% (1.9% - 2.0%)	2,672	0.8% (0.6% - 1.0%)	1.3% (1.3% - 1.4%)	1.4% (1.3% - 1.4%)
7	1,412	1.2% (1.0% - 1.6%)	2.8% (2.7% - 2.8%)	2.3% (2.3% - 2.4%)	1,315	1.0% (0.7% - 1.2%)	1.5% (1.5% - 1.5%)	1.5% (1.5% - 1.6%)
8	677	1.5% (1.1% - 1.9%)	3.3% (3.2% - 3.3%)	2.7% (2.7% - 2.8%)	591	1.1% (0.8% - 1.5%)	1.7% (1.7% - 1.7%)	1.8% (1.7% - 1.8%)
9	226	1.7% (1.2% - 2.6%)	3.8% (3.8% - 3.9%)	3.2% (3.1% - 3.2%)	177	1.1% (0.8% - 1.8%)	1.9% (1.8% - 2.0%)	2.0% (2.0% - 2.0%)
10		( - )	4.4% (4.3% - 4.5%)	3.6% (3.6% - 3.7%)		( - )	2.1% (2.1% - 2.2%)	2.3% (2.2% - 2.3%)
11		( - )	4.9% (4.8% - 5.0%)	4.1% (4.0% - 4.2%)		( - )	2.4% (2.4% - 2.5%)	2.6% (2.6% - 2.7%)
12		( - )	5.5% (5.4% - 5.6%)	4.6% (4.5% - 4.7%)		( - )	2.7% (2.6% - 2.8%)	3.0% (2.9% - 3.0%)
13		( - )	6.0% (5.9% - 6.1%)	5.1% (5.0% - 5.2%)		( - )	3.1% (3.0% - 3.2%)	3.4% (3.3% - 3.5%)
14		( - )	6.6% (6.4% - 6.8%)	5.6% (5.5% - 5.7%)		( - )	3.5% (3.4% - 3.7%)	3.8% (3.7% - 3.9%)
15		( - )	7.3% (7.0% - 7.6%)	6.2% (6.0% - 6.5%)		( - )	( - )	4.5% (4.2% - 4.7%)

Cumulative revision rate with 95% confidence intervals

Rate is only reported for times where more than 40 remain at risk

## Patient Reported Outcomes

### PROMs Analysis

Comprising PROMs data up to and including: 31/12/2016

Oxford Hip Score (0 - 48)	R3 Cementless Cup	All NJR Cementless Cups	All THR in NJR
Paired Records	6,006	156,545	247,104
PreOp score	18.6 (18.4 - 18.8)	18.5 (18.5 - 18.5)	18.1 (18.1 - 18.1)
6 month score	40.4 (40.2 - 40.6)	39.8 (39.8 - 39.9)	39.2 (39.2 - 39.3)
6 month score (adjusted)	40.0 (39.8 - 40.2)	39.5 (39.5 - 39.5)	39.2 (39.2 - 39.3)
Health gain	21.8 (21.6 - 22.1)	21.3 (21.3 - 21.4)	21.1 (21.1 - 21.2)
Health gain (adjusted)	21.9 (21.7 - 22.1)	21.4 (21.4 - 21.4)	21.1 (21.1 - 21.2)
p value (adjusted health gain)	< 0.001		
Score Improved	97.2%	97.2%	97.1%

EQ-5D Index (-0.59 - 1.00)	R3 Cementless Cup	All NJR Cementless Cups	All THR in NJR
Paired Records	5,569	144,340	226,875
PreOp score	0.373 (0.365 - 0.382)	0.371 (0.369 - 0.372)	0.357 (0.356 - 0.359)
6 month score	0.810 (0.804 - 0.817)	0.801 (0.800 - 0.803)	0.789 (0.788 - 0.790)
6 month score (adjusted)	0.805 (0.799 - 0.811)	0.795 (0.794 - 0.796)	0.789 (0.788 - 0.790)
Health gain	0.437 (0.428 - 0.446)	0.431 (0.429 - 0.433)	0.432 (0.431 - 0.433)
Health gain (adjusted)	0.448 (0.442 - 0.454)	0.438 (0.437 - 0.439)	0.432 (0.431 - 0.433)
p value (adjusted health gain)	< 0.001		
Score Improved	89.6%	89.5%	89.2%

EQ-VAS (0 - 100)	R3 Cementless Cup	All NJR Cementless Cups	All THR in NJR
Paired Records	5,318	138,951	217,695
PreOp score	65.7 (65.1 - 66.3)	65.7 (65.6 - 65.8)	65.3 (65.2 - 65.4)
6 month score	78.7 (78.2 - 79.1)	77.5 (77.4 - 77.6)	76.6 (76.6 - 76.7)
6 month score (adjusted)	78.1 (77.7 - 78.5)	76.9 (76.8 - 77.0)	76.6 (76.6 - 76.7)
Health gain	13.0 (12.4 - 13.6)	11.8 (11.7 - 11.9)	11.3 (11.2 - 11.4)
Health gain (adjusted)	12.8 (12.4 - 13.2)	11.6 (11.5 - 11.7)	11.3 (11.3 - 11.4)
p value (adjusted health gain)	< 0.001		
Score Improved	68.9%	66.1%	65.0%

Adjusted scores correspond to the NHS Digital version 3 casemix-adjustment model. This includes patient variables from PROMs, HES, and IMD data.

Responses are not case mix adjusted

**Patient Reported Outcomes**



## Disclaimer

### Disclaimer

The National Joint Registry (NJR) produces this report through its contractor, Northgate Public Services (UK) Limited (NPS) using data collected, collated and provided by third parties. As a result, neither NJR or NPS takes any responsibility for the accuracy, currency, reliability and correctness of any data used or referred to in this report, and disclaims all warranties in relation to such data, in each case to the maximum extent permitted by legislation.

To the maximum extent permitted by law, neither the NJR or NPS shall have any liability (including but not limited to liability by reason of negligence) for any loss, damage, cost or expense incurred or arising by reason of any person using or relying on the data within this report and whether caused by reason of any error, omission or misrepresentation in the report or otherwise. This report is not to be taken as advice. Third parties using or relying on the data in this report do so at their own risk and will be responsible for making their own assessment and should verify all relevant representations, statements and information with their own professional advisers.

## Appendix A: Component List

Catalogue Number	Description	Number Implanted
71330004	R3 HEMI 3 HOLE ACET SHEL MM52	3
71330005	R3 HEMI 3 HOLE ACET SHEL MM54	2
71330006	R3 HEMI 3 HOLE ACET SHEL MM56	2
71330007	R3 HEMI 3 HOLE ACET SHEL MM58	2
71331854	R3 0 HOLE ACET SHELL 54MM	1
71331940	R3 3 HOLE HA CTD ACET SHELL 40MM	5
71331942	R3 3 HOLE HA CTD ACET SHELL 42MM	3
71331944	R3 3 HOLE HA CTD ACET SHELL 44MM	9
71331946	R3 3 HOLE HA CTD ACET SHELL 46MM	352
71331948	R3 3 HOLE HA CTD ACET SHELL 48MM	1,951
71331950	R3 3 HOLE HA CTD ACET SHELL 50MM	3,462
71331952	R3 3 HOLE HA CTD ACET SHELL 52MM	4,657
71331954	R3 3 HOLE HA CTD ACET SHELL 54MM	3,460
71331956	R3 3 HOLE HA CTD ACET SHELL 56MM	2,618
71331958	R3 3 HOLE HA CTD ACET SHELL 58MM	1,493
71331960	R3 3 HOLE HA CTD ACET SHELL 60MM	711
71331962	R3 3 HOLE HA CTD ACET SHELL 62MM	262
71331964	R3 3 HOLE HA CTD ACET SHELL 64MM	104
71331966	R3 3 HOLE HA CTD ACET SHELL 66MM	11
71331968	R3 3 HOLE HA CTD ACET SHELL 68MM	1
71332246	R3 0 HOLE HA CTD ACET SHELL 46MM	34
71332248	R3 0 HOLE HA CTD ACET SHELL 48MM	285
71332250	R3 0 HOLE HA CTD ACET SHELL 50MM	505
71332252	R3 0 HOLE HA CTD ACET SHELL 52MM	703
71332254	R3 0 HOLE HA CTD ACET SHELL 54MM	466
71332256	R3 0 HOLE HA CTD ACET SHELL 56MM	421
71332258	R3 0 HOLE HA CTD ACET SHELL 58MM	210
71332260	R3 0 HOLE HA CTD ACET SHELL 60MM	98
71332262	R3 0 HOLE HA CTD ACET SHELL 62MM	25
71332264	R3 0 HOLE HA CTD ACET SHELL 64MM	3
71335548	R3 3 HOLE ACET SHELL MM48	2
71335550	R3 3 HOLE ACET SHELL 50MM	6
71335552	R3 3 HOLE ACET SHELL 52MM	1
71335554	R3 3 HOLE ACET SHELL 54MM	1
71335556	R3 3 HOLE ACET SHELL 56MM	1
71338663	R3 MULTI HOLE ACETABULAR SHELL 48MM	6
71338664	R3 MULTI HOLE ACETABULAR SHELL 50MM	9
71338665	R3 MULTI HOLE ACETABULAR SHELL 52MM	14
71338666	R3 MULTI HOLE ACETABULAR SHELL 54MM	18
71338667	R3 MULTI HOLE ACETABULAR SHELL 56MM	8
71338668	R3 MULTI HOLE ACETABULAR SHELL 58MM	8
71338669	R3 MULTI HOLE ACETABULAR SHELL 60MM	2
71338671	R3 MULTI HOLE ACETABULAR SHELL 62MM	3
71338733	R3 HA MULTI-HOLE ACET SHELL 48MM	7
71338734	R3 HA MULTI-HOLE ACET SHELL 50MM	17
71338735	R3 HA MULTI-HOLE ACET SHELL 52MM	20

## Appendix A: Component List

Catalogue Number	Description	Number Implanted
71338736	R3 HA MULTI-HOLE ACET SHELL 54MM	10
71338737	R3 HA MULTI-HOLE ACET SHELL 56MM	12
71338738	R3 HA MULTI-HOLE ACET SHELL 58MM	11
71338739	R3 HA MULTI-HOLE ACET SHELL 60MM	2
71338740	R3 HA MULTI-HOLE ACET SHELL 62MM	3
71338741	R3 HA MULTI-HOLE ACET SHELL 64MM	3
71338742	R3 HA MULTI-HOLE ACET SHELL 66MM	1
71338743	R3 HA MULTI-HOLE ACET SHELL 68MM	1
71338745	R3 HA MULTI-HOLE ACET SHELL 72MM	1
71339452	R3 MAX HEMI 3 HOLE MM52	4