## ACUFEX® Surgical Instruments

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### Intended Use

ACUFEX Surgical Instruments are intended for use in orthopedic and general surgery.

### ⚠️ Warnings

- This product is shipped non-sterile. It must be disassembled, cleaned, and sterilized before the first use. It must be disassembled, cleaned, and sterilized before every subsequent use.**

- It is the surgeon's responsibility to be familiar with the appropriate surgical techniques prior to use of this device.**

- Read these instructions completely prior to use.**

### Precautions

The U.S. Federal law restricts this device to sale by or on the order of a physician.

- Prior to use, inspect the device to ensure it is not damaged. Do not use a damaged device.
- As with any surgical instrument, careful attention should be exercised to ensure that excessive force is not placed on the instrument. Excessive force can result in instrument failure.
- Do not use these instruments as levers for manipulating hard tissue or bone. Excessive force should not be applied to the instrument when manipulating soft tissue, bone, or hard objects. Misuse of these instruments may result in bent distal tips or jaws, and dull or uneven cutting edges.
- Careful attention must be paid to asepsis and avoidance of anatomical hazards.

### Cleaning

**Notes:**

- Pay careful attention to cleaning devices with challenging design features. Challenging design features can include, but not limited to, suction levers, stopcocks, interfaces, cannulations, holes, blind holes, crevices, hinges, mating surfaces, etc.
- Low-sudsing, neutral 6.0–8.0 pH, enzymatic detergents are recommended.
- Do not use detergents above 11.0 pH.
- Use deionized water for washing and rinsing.

#### Manual Cleaning

- Disassemble the instrument into its component parts.
- Ensure that suction levers or stopcocks are in the full open position or disassembled.
- Immerse the instrument and soak for a minimum of five minutes in enzymatic detergent and deionized water.
- Remove soil from challenging design features with cleaning brushes. Scrub interfaces, cannulations, and holes with a light fitting brush using a twisting action. If possible, retract or move components to access and clean these areas. Scrub crevices and around hinged/mating surfaces with a brush.
- Sonicate the instrument in warm enzymatic detergent and deionized water in its fully open position for a minimum of 15 minutes.
- Rinse thoroughly with warm deionized water, making sure to irrigate the challenging design features. Blind holes should be repeatedly filled and emptied.
- Check the instrument for any remaining soil. Repeat cleaning if soil is visible.

English

## Automatic Cleaning

- Disassemble the instrument into its component parts.
- Ensure that suction levers or stopcocks are in the full open position or disassembled.
- Immerse in enzymatic detergent and deionized water. Remove soil from challenging design features with cleaning brushes. Scrub interfaces, cannulations, and holes with a light fitting brush using a twisting action. If possible, retract or move components to access and clean these areas. Scrub crevices and around hinged/mating surfaces with a brush.
- Sonicate in fully open position for a minimum of 10 minutes in warm enzymatic detergent and deionized water.
- Rinse with warm deionized water.

- Put the devices through an automatic washing cycle. Follow the instructions included with the automatic washer for instrument load and configuration.

**Minimum Parameters:**

- Use deionized water
- 5 minute cold prewash
- 5 minute enzymatic wash at 109° F (43° C)
- 5 minute detergent wash at 131° F (55° C)
- 1 minute rinse at 113° F (45° C)

- Thermal disinfection 196° F (91° C) for 1 minute minimum.

### Sterilization

**Note:** If a Smith & Nephew instrument tray is used for sterilization, refer to the instructions for use provided with the instrument tray for the recommended sterilization parameters.

Endoscopic instruments may be sterilized by ethylene oxide or by steam. Time and temperature parameters required for sterilization vary according to type of sterilizer, cycle design, and packaging material. Please review the sterilizer manufacturer's sterilization instructions, or hospital procedures, prior to sterilization.

#### Steam Sterilization

The minimum recommended parameters for steam sterilization are as follows:

**Standard Gravity Steam Cycle**

Temperature: 250° F (121° C)

Minimum exposure time: 30 minutes

**High Temperature Gravity Steam Cycle**

Temperature: 270° F (132° C)

Minimum exposure time: 10 minutes

**Prevacuum Cycle**

Temperature: 270° F (132° C)

Minimum exposure time: 3 minutes

**World Health Organization (WHO) Steam Cycle**

Temperature: 273–280° F (134–138° C)

Minimum exposure time: 18 minutes

#### Ethylene Oxide Sterilization

The minimum recommended parameters for ethylene oxide sterilization are as follows:

Temperature: 122–140° F (50–60° C)  
Humidity: 60–100%  
EtO concentration: 580–730 mg/l  
Exposure time: 4 hours  
Aeration time: 12 hours

**Note:** Sterilizing in liquid solutions is not recommended.

### Warranty

This product is warranted to be free from defects in material and workmanship.

### For Further Information

If further information on this product is needed, please contact Smith & Nephew Customer Service at +1-800-343-5717 in the U.S., or your authorized representative.

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